



**Glenammer**  
Laboratory Test Sieves

# Gem Series

## Advanced Shaker



## Features

Users can adjust the amplitude from 0-9 steps in order to control the sieving power and shaking intensity.

Introduces the function of an interval operation. User can either run the shaker continuously or with intervals.

GEM Advance sieve shaker allows users to retain 5 sets of memory for restore. Users can simply select the pre-named number and bring recall all the settings.

# Gem Advanced

### Capacity

8 x 200mm Diam (100mm, 150mm, 3" and 8") sieves plus lid and receiver.  
Extended rods available.

### Net Weight

46Kg

### Max sample weight

3Kg

### Sieve Motion

3D Electromagnetic

### Amplitude

0-3mm (max amplitude depends on loading). Set digitally using 0-9 steps

### Electrical Supply

230 Volt, 1 ph, 50 Hz.  
Other voltages and frequencies available

### Interval

On / Off interval operation

### Timer

Digital: 0sec - 99mins 59secs

### Net W X H X D

400mm x 200mm x 410mm

### Memory function

Store and recal up to 5 set operaton settings