



Glenammer
Laboratory Test Sieves

GEM 200

Analogue & Digital Shakers



Features

GEM 200 Digital has a digital display panel with an installed digital timer, users can set up the desired operation sieving time between 0 sec to 99 mins 59 secs.

The system features a pause/start option, which enables users to operate the shaker with intervals.

GEM 200 Analogue has an analogue timer to adjust the operating time between 0min-60min, users can simply turn the dial to the desired time to start the operation.



GEM 200 Analogue & Digital



Capacity

8 X 200mm Dia. (100mm, 150mm, 3" and 8") sieves plus lid and receiver.



Net Weight

20Kg



Max sample weight

3Kg



Sieve Motion

Electromagnetic



Speed

3,000/min at 50 Hz



Electrical Supply

230 Volt, 1 ph, 50 Hz, Input power 0.045 kw, current 0.20 amps



Height

140mm



Timer

Digital: 0sec - 99mins 59secs
Analogue: 0 - 60min



Net W x H x D

240mm x 140mm x 240mm



Type

8 x 200mm dia
(100mm, 150mm, 8")
sieves plus lid and receiver